

Safersonic Plus

**Sterile Transducer Cover with Adhesive
+ Safergel Ultrasound Transmission Gel**



Couples advantages of Safersonic Transducer Cover with Adhesive:

- *No Gel at the Transducer*
- *No Artifacts*

and Safergel Ultrasound Transmission Gel:

- *NO PRESERVATIVES*
- *Precise and easy application*



unites advantages of **Safersonic** **Transducer Cover with Adhesive:**

- No Rubber Bands
- No Artifacts
- No Gel at the Transducer
- Latexfree

and **Safergel** **Ultrasound Transmission Gel:**

- NO PRESERVATIVES
- Bubble Free Delivery
- Optimized Viscosity
- Easy, precise application
- „TWIST`N USE“ SEAL

The composition of the Innovations-award winning tube (ETMA Award 2015) allows EO-sterilization of cover together with Gel!

For IFU and videos please visit: www.safersonic.com

ProdNbr.	Item	Dim.	Gel	Qty
SAF.612.18030.WG	Safersonic Single Plus	18 x 30 cm	15g Safergel	50
SAF.612.18050.WG*	Safersonic Medium Plus	18 x 50 cm	15g Safergel	50
SAF.612.18120.WG	Safersonic Conti Plus	18 x 120 cm	15g Safergel	50
SAF.612.18250.WG	Safersonic Long Plus	18 x 250 cm	15g Safergel	40

* on request

DISTRIBUTOR

Safersonic Medizinprodukte Handels-GmbH

Herrengasse 5 | A - 3370 Ybbs an der Donau | Tel.: +43 (0) 7412 55772 | E-mail: info@safersonic.com